



umano medical

## AN ALL-IN-ONE HOSPITAL BED FOR MULTIPLE CARE ENVIRONMENTS

med-surg | bariatric care | palliative care



 **ook snow** » ALL + **Pilot Drive**™

One-step powered Infinite™ width expansion system, increased weight capacity, and low height all come together to improve staff efficiency, patient comfort, and safety.

Featuring Pilot Drive™ – the first powered drive system to offer an unmatched combination of low height, underbed clearance, and manoeuvrability.

### VERSATILITY • LOW HEIGHT • MOBILITY



**MEDICAL-SURGICAL**  
Low height, ease of use, and ease of cleaning makes this bed ideal for med-surg environments.



**BARIATRIC CARE**  
1060 lb weight capacity, width expansion, and an appealing design encourages patient autonomy and dignity for bariatric care.



**PALLIATIVE CARE**  
Also known as the Cuddle Bed, it is the only bed on the market certified to be used by two people.



**PILOT DRIVE**  
Added mobility combined with enhanced versatility makes the ook snow ALL with Pilot Drive and integrated O2 bottle holders a unique solution for your bariatric patients in intensive care.



STANDARD FEATURES

- 11" low height with 5" dual casters
Infinite™ powered width adjustment
Integrated scale system
Smart screen, nurse controls and patient control lockout on footboard
Integrated bed exit detection system
3 levels of sensitivity
Bed exit side view lights
Smart Nightlights™
Remote bed exit alarm connectivity
Boostless™ Backrest System
One hand release siderails
Liftable and flat deck covers
Battery back-up
Compatibility with mattresses within 48" width
Programmable lowest height
Nurse and patient controls integrated into siderails
Integrated bed extender (35" x 4")
CPR release handle
Trendelenburg angle displayed on footboard
Backrest angle indicator
Patient restraint loops
Integrated pump holder
120-volt / 5A insulated outlets (at foot end)
Communication port (for bed software updates)
Programmable Auto-Contour mode
One-touch button groups
30-degree head of bed
One-button chair position
Flat bed position
Trendelenburg / Reverse Trendelenburg
Vascular foot position

OPTIONAL FEATURES

- 12.7" low height 6" dual casters
Pilot Drive™ with or without O2 bottle holders (compatible with size E bottles)
Downward motion interrupt
Bed status monitoring
Optional patient surface length (80" or 84")
120-volt/5A insulated outlets with USB mobile device power outlets
Higher siderails (for 8" mattress)
Foot end back-up control on the frame
Roller type wall protection bumpers
Additional nurse call port
One-touch button groups
Backlit patient controls
Nurse call
Different styling options available

ACCESSORIES

- Patient helper system (300 lb capacity)
Upright oxygen bottle holder
IV poles (permanent & removable)
4" bed extender cushion
Removable patient controls
Transport cart

SPECIFICATIONS



Table with 2 columns: Specification Name and Value. Includes rows for Product model (200-0000), Length (Overall length 80"/84", Extended length with Pilot Drive), Overall width (Siderails down/up), Safe Working Load (1060 lb / 481 kg), Height range (to deck top), Patient deck positioning (Backrest, Knee gatch, Trend./Reverse Trend., Boostless™), Patient deck dimensions (Width, Length, Extended Length), Caster diameter (with Pilot Drive).

Compliant with standards IEC 60601-1 and IEC 60601-2-52. Certified for use by two people (referring to one patient and one bystander.)

ELECTRONICS

Table with 2 columns: Electronics Specification and Value. Includes rows for Volts (120 V AC), Ampere rating (4.25A), Current leakage (Compliant with IEC 60601-1), Frequency (60 Hz), Hospital grade plug (Yes), Battery type (Sealed rechargeable lead acid 1.2Ah), Battery type (Pilot Drive) (Sealed rechargeable lead acid 16.5Ah), Scale system accuracy (1% with a maximum difference of 2.2 lb / 1 kg), Scale system accuracy (EC) (OIML Class III ± 1 kg).

Umano Medical reserves the right to change specifications without notice.

Umano Medical Inc.

230, boulevard Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll free number : 1 844-409-4030
information@umanomedical.com
umanomedical.com

