



## Breezy Relax Transit Wheelchair

PRODUCT CODE NB0530 | WHEELCHAIRS

### Features

- Gas assisted, tilt-in-space and recline functions provide optimum postural positioning
- Removable adjustable height arm feature padded side panels and adjustable armpad length
- Removable swingaway elevating legrests and calf pads
- Height adjustable padded contoured backrest
- Depth adjustable padded seat with inbuiltommel
- Removable padded headrest with adjustable side supports
- Tilt angle -2° to 20°
- Recline angle 0° to 40°
- Drum brake option allows attendant to control brakes from the push handles
- Alloy construction with water resistant upholstery
- 16" rear wheels with solid tyres

### Additional Information

- ARTG Number - 199075

### Specifications

- Seat Width: 410mm
- Seat Depth: 480-500mm
- Seat Height: 410-520mm
- Back Height: 500-650mm
- Arm Height: 180-290mm
- Push Handle Height: 980-1160mm
- Width Between Arms: 360mm
- Overall Width: 640mm
- Gross Weight: 38kg
- Net Weight: 22kg
- Max User Weight: 140kg