

# Alpha Lifecare



## Avalo<sup>®</sup> Medical Carts

# Setting the Standard for Efficient Workflow and Reliable Performance

Avalo® Medical Carts careful design, functionality, and logical organisation, enhance user efficiency and support the delivery of superior patient care. Its versatility makes it ideal for various care applications, such as hospitals, surgery centers, physicians' offices, assisted living facilities, and more.



## CONVENIENT BUNDLES

Experience the ease of purchasing with our thoughtfully curated bundles that combine a cart and accessory package for various care applications.



## BESPOKE SOLUTIONS

With fully configurable options, you have the freedom to design your own cart, customising drawer configurations and accessories to seamlessly support your facility's unique workflow and storage requirements.



## FEATURES

- » Precision-welded steel frame with durable, high-impact panels
- » Lightweight design and profile simplify movement and steering
- » Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- » Innovative locking systems including key, keyless or quick-access breakaway locking
- » Available in Standard, Intermediate, or Compact height
- » Selection of three drawer depths, 3", 6" and 10"
- » Wide selection of accessories and organisation systems
- » Expansive work area with optional slide-out work surface
- » Delivered fully assembled, no tools required
- » 10-year limited warranty – structure
- » 2-year limited warranty - key locking mechanism, castors, draw runner system.

## SPECIFICATIONS

Cart	Height (cm)	Depth (cm)	Width (cm)
Standard	110	61	79
Intermediate	100	61	79
Compact	91	61	79
Drawer			
3 Inch	7	45	55
6 Inch	16	45	55
10 Inch	25	45	55

# SECURITY FEATURES

## & LOCK SYSTEM OPTIONS

Avalo carts stand out as the premier choice in medical cart solutions, offering an unparalleled advantage by providing a comprehensive array of security features for customisation.



## LOCK SYSTEM CHOICES



### AutoLock™

Automatically secures a motionless unlocked medical cart with a user-defined relock time delay from 1-99 minutes. AutoLock permits up to 5,000 individual user codes to further enhance cart access tracking and security. SLA battery operated. **AutoLock security upgrades available.**



### Keyless Access with Manual Relocking

Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medical cart, simply push in the lock lever. Clear visual indication of lock status is apparent from a distance.



### Key Access with Manual Relocking

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides your medical cart advanced security and the flexibility to change keys simply and easily.



### Breakaway Lock

Avalo Emergency Cart's unique breakaway locking handle is standard on every model, offers the security you require and ensures quick access to code cart supplies.

## AUTOLOCK™ SECURITY UPGRADES



### Standard Host Software

Using Standard Host Software authorised users can add, change or delete an access code for specific users.



### Proximity Access

An advanced proximity card reader can be selected as an additional security access requirement.



### Keyless Narcotics Storage

The optional keyless narcotics storage feature greatly augments security to further optimise control and accountability.

**Available for Avalo AC Cart only**

# EMERGENCY CART

Designed for reliable performance and improved workflow in emergency departments, the Avalo® Emergency Cart boasts a unique breakaway lock for optimal security and easy access. With stable manoeuvrability and accessories, it enhances user efficiency and supports code team performance.

**WELL SUITED FOR:**

Acute Care Hospitals, Extended Care facilities and Nursing Homes, Residential Care Facilities Rehab and Psych Environments, Prompt Care Clinics, Dental Offices, Physician Offices.

# PAEDIATRIC CRASH CART

Designed for use with the proven clinical methodology for treating life-threatening emergencies in paediatric patients. With colour-coded drawers, doctors and clinicians can swiftly identify critical supplies when every second matters. Contents are secured with a breakaway locking, providing immediate access in urgent situations with a simple swift motion.

**WELL SUITED FOR:**

Acute Care Hospitals, Children's Hospitals, Paediatric Emergency Departments.

## EMERGENCY CART BUNDLE

**STEP 1: Choose a Base Cart Option**



**1. STANDARD HEIGHT**

- Includes:
- » Breakaway lock
  - » 3 x 3" drawer
  - » 2 x 6" drawer
  - » 1 x 10" drawer
  - » 2 x A&E mounts
  - » Emergency label red

Red | GC3700



**2. INTERMEDIATE HEIGHT**

- Includes:
- » Breakaway lock
  - » 2 x 3" drawer
  - » 2 x 6" drawer
  - » 1 x 10" drawer
  - » 2 x A&E mounts
  - » Emergency label red

Red | GC3600



**3. COMPACT HEIGHT**

- Includes:
- » Breakaway lock
  - » 3 x 3" drawer
  - » 1 x 6" drawer
  - » 1 x 10" drawer
  - » 2 x A&E mounts
  - » Emergency label red

Red | GC3500

**STEP 2: Add Optional Accessory Package | GZ4540**



Tamper Evident Seals (100)



Oxygen Tank Holder



Defibrillator Shelf



IV Pole



Cardiac Board



3" Drawer Dividers (Set of 2)

## PAEDIATRIC BUNDLE

**STEP 1: Choose the Base Cart Option**

**STANDARD HEIGHT**

- Includes:
- » Breakaway locking system
  - » 8 x 3" drawer
  - » 1 x 6" drawer with dividers
  - » 2 x A&E mounts

Light crème | GC3800



**STEP 2: Add Optional Accessory Package | GZ4540**



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



Defibrillator Shelf



IV Pole



3" Drawer Dividers (Set of 2)



# ANAESTHESIA CART

Designed to meet the unique needs of your anaesthesia department, bringing together function and technology in perfect harmony. Choose from (3) lock systems, including options for secondary proximity access for added security. Flexible divider systems and removable drawer trays allow for smart organisation of medications and all anaesthesia-related supplies.

**WELL SUITED FOR:**

Acute Care Hospitals, Skilled Nursing Facilities, Outpatient Surgery Centers, Dental Surgery Office.

# ISOLATION CART

Designed with storage space a priority, this cart allows you to keep all your isolation supplies neatly organised and easily accessible with drawers that open fully so you can reach contents in the rear effortlessly. Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention. Drawers stay securely closed during motion.

**WELL SUITED FOR:**

Extended Care Facilities and Nursing Homes, Acute Care Hospitals, Rehab and Psych Environment.

## ANAESTHESIA BUNDLE

**STEP 1: Choose a Base Cart Option**



**1. STANDARD HEIGHT**

- Includes:
- » Keyless locking system
  - » 3 x 3" drawer
  - » 2 x 6" drawer with dividers
  - » 1 x 10" drawer with dividers
  - » 2 x A&E mounts
  - » Anaesthesia cart label, extreme green

Extreme green | GC4200



**2. INTERMEDIATE HEIGHT**

- Includes:
- » Keyless locking system
  - » 2 x 3" drawer
  - » 2 x 6" drawer with dividers
  - » 1 x 10" drawer with dividers
  - » 2 x A&E mounts
  - » Anaesthesia cart label, extreme green

Extreme green | GC4100



**3. COMPACT HEIGHT**

- Includes:
- » Keyless locking system
  - » 3 x 3" drawer
  - » 1 x 6" drawer with dividers
  - » 1 x 10" drawer with dividers
  - » 2 x A&E mounts
  - » Anaesthesia cart label, extreme green

Extreme green | GC4000

**STEP 2: Add Optional Accessory Package | GZ4550**



Writing Surface



Waste with Lid



3" Drawer Dividers (Set of 3)



Storage Tower



Tilt 4 Bin Unit



Tilt Bin 5 Bin Unit

## ISOLATION BUNDLE

**STEP 1: Choose the Base Cart Option**

**INTERMEDIATE HEIGHT**

- Includes:
- » No locking system
  - » 3 x 10" drawer with dividers
  - » 2 x A&E mounts
  - » Isolation cart label


Extreme yellow | GC4300



**STEP 2: Add Optional Accessory Package | GZ4560**



3" Drawer Dividers (1)



Writing Surface



Waste Bin with Lid



# PROCEDURE CART

Designed with ample storage capacity this cart is adaptable to various healthcare applications including trauma, airway, dressing, cast, or a general supply cart. Security is prioritised with three (3) proven lock systems, including an advanced keyless entry with auto-relock feature.

**WELL SUITED FOR:**

Extended Care Facilities and Nursing Homes as a Medication Cart, Acute Care Hospitals, Residential Care Facilities, Rehab and Psych Environments  
Prompt Care Clinics, Dental Offices, Physician Offices.

# MEDICATION CART

The perfect solution for hospital environments, combining cutting-edge design with ample capacity and easy manoeuvrability. Customise your medication control with a choice of cassette size, ensuring highly secure and organised medication management. Despite its compact design, these carts offer generous storage and remain lightweight, making it easy for healthcare personnel to handle.

**WELL SUITED FOR:**

Nursing Home/Extended Care Facility, Assisted Living, Continuing Care Community, Hospital.

## PROCEDURE BUNDLE

**STEP 1: Choose a Base Cart Option**



**1. INTERMEDIATE HEIGHT**

- Includes:
- » Choice of Key Lock or Keyless
  - » 2 x 3" drawer
  - » 2 x 6" drawer
  - » 1 x 10" drawer
  - » 2 x A&E mounts
  - » Dark creme procedure cart label red

Dark crème with key lock | GC4600  
Dark crème with keyless | GC4700



**2. COMPACT HEIGHT**

- Includes:
- » Choice of Key Lock or Keyless
  - » 1 x 3" drawer
  - » 2 x 6" drawer
  - » 1 x 10" drawer
  - » 2 x A&E mounts
  - » Dark creme procedure cart label red

Dark crème with key lock | GC4400  
Dark crème with keyless | GC4500

**STEP 2: Add Optional Accessory Package | GZ4560**



Waste Bin with Lid



3" Drawer Divider (1)



Writing Surface

## YOUR CHOICE OF 2 MODELS



**ACM**

- » Keyless locking system
- » 2 x 3" drawer
- » 1 x 10" drawer
- » 1 x 5.5" two tier cassette package
- » 1 x 5" three tier cassette package
- » 2 x 3" utility drawer
- » 1 x 6" utility drawer
- » 2 x 10" utility drawer

Dark crème | GC5100



**ACS9**

- » Keyless locking system
- » 2 x 3" main drawer
- » 1 x 10" main drawer
- » 2 x 8" ac two tier cassette package

Dark crème | GC5000



# OPTIONS & ACCESSORIES\*

Purchasing the accessory packages in bundles equips you with a well-rounded set of tools and features. Should you require additional items, you have the option to add to more accessories your cart.



Cardiac Board



IV Pole



Defibrillator Shelf



Oxygen Tank Holder



Divider Drawer Tray



Drawer Dividers



Bulk Storage Bins



Tilt Bin Organisers



Side Organiser



Waste Containe



Sharps/Gloves Holder



Storage Container



Internal Lock Box



Writing Surface



Drawer Tag Holder



Tracking Castor



Accessory Bridge



Clear Top Mat



Plug Outlet Strip  
UL1363A Rated



Computer Arm

\*Extended lead times and additional charges may apply.

# BUILD YOUR OWN CART


Custom configure an Avalo Medical Cart to meet your clinical workflow requirements.

## STEP 1

Design your cart online:  
[configurator.capsahealthcare.com/avalo](https://configurator.capsahealthcare.com/avalo)

## STEP 2

Save your customised configuration results as a PDF and send to [info@alphalifecare.com.au](mailto:info@alphalifecare.com.au) for a quote.



\*Extended lead times and additional charges may apply.

## Caring Further

Level 45  
88 Walker St  
North Sydney, NSW 2060

1300 930 930  
[info@alphalifecare.com.au](mailto:info@alphalifecare.com.au)  
[alphalifecare.com.au](https://alphalifecare.com.au)

### Australia wide locations

Sydney | Newcastle | Brisbane | Melbourne  
Canberra | Griffith | Tasmania

© Copyright Alpha Lifecare Pty Ltd 2024, all rights reserved.

Alpha Lifecare

