



## Introducing the Casero bed series

*Meet the bed designed for  
your needs and your home*



*Delivering home style luxury beds  
across all care settings*







## Features of the future

The Casero bed series showcases sophisticated, stylish and intelligent solutions to suit both client's needs and their homes.

Developed especially for homecare and aged care environments, the Casero series offers 4 different bed styles covering aged care and community clients, alongside additional and highly functional features to reduce falls risks as well as options for bariatric users and couples.

The Casero series is modular and configurable to your client's requirements and includes classic and contemporary designs, upholstery choices and high quality features.

The European designed 4-part telescopic side rail option offers versatility, support, safety, protection – and overall peace of mind. With an extensive range of options and accessories, the Casero series truly reimagines and merges perfectly the functional needs of our clients and the style that is important to them.













## Casero styles & finishes

The Casero range is modular and configurable to your client's requirements and includes classic and contemporary designs, upholstery choices and high quality features.

Headboard styles		Colour options	
The Oxford		Oatmeal	
		Granite	
The Belgrave		Slate	
		Ebony (vinyl)	

\*All coverings available in each headboard style. Colours may vary due to printer being used.



# Casero features & accessories

## The Duo

Width and length adjustable  
Height range 184mm - 765mm  
Max user weight 250kg  
Single (BD6000)  
King Single (BD6100)

---

## The Mode

Disassembles for easy transportation ideal for homes with restricted accessibility  
Low height range ideal for falls risk clients  
Height range 110mm - 740mm  
Max user weight 200kg  
Single (BD6200)  
King Single (BD6300)

---

## The Vogue

Low height range ideal for falls risk clients  
Height range 110mm - 740mm  
Max user weight 225kg  
Single (BD6400)  
King Single (BD5000)

---

## The Plus

Width and length adjustable  
Suitable for couples or clients requiring greater sleeping area or increased weight capacity  
Height range 250mm - 660mm  
Max user weight 450kg  
Double (BD6600)  
Queen (BD6700)

## Features

Electric Hi/Low Functions  
Electric Back Elevation  
Electric Knee Break  
Electric Trendelenburg\*  
Built-in 200mm length extension

## Accessories

Self-help pole (BZ0325)  
3/4 length horizontal side rails (BZ0443)  
Telescopic 4 part side rail (BZ0302)  
Battery back up system (BZ0335)\*  
Nurse control panel (BZ0337)\*  
Under bed light (BZ0336)\*

## Warranties

Steel frame - 10 years  
Electrics - 5 years  
Hand control - 12 months  
Upholstered surround incl head and foot boards - 12 months

\* not available on The Plus model



# Making lives better

The Casero bed series showcases sophisticated, stylish and intelligent solutions to suit both client's needs and their homes

---

## Australia wide locations

Sydney | Newcastle | Brisbane  
Melbourne | Canberra | Griffith

1300 930 930  
info@alphalifecare.com.au  
**alphalifecare.com.au**

GPO Box 3390  
Sydney NSW 2001  
Australia

---