



## Juvo Mobile, Swingaway Legrests - 460mm Width

PRODUCT CODE AE0510 | MOBILE SHOWER COMMODES

### Features

- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet for toileting
- Padded commode seat (available as open or closed front)
- Swingback arms with padded armpads for side transfers
- Armrests can be adjusted to two different heights
- Adjustable safety arms for additional security
- PVC slatted backrest
- Adjustable height swingaway legrests feature a unique unobtrusive mechanism design to prevent skin tears
- Mounted on 4 × 125mm locking stainless steel castors
- Durable stainless steel frame
- Tested and conforms to Australian Standards AS/NZS 3973: 2009

### Also Available

- AE1050 - Attendant Propelled and Height Adjustable

### Sold Separately

- Padded commode seats (open or closed front) are additional

### Additional Information

- ARTG Number - 191619

### Specifications

- Seat Width - 460mm
- Seat Depth - 430mm
- Seat Height - 565mm
- Armrest Height - 200-240mm
- Backrest Height - 340mm
- Overall Width - 540mm
- Max User Weight - 200kg