



MOVI® Mattress Solutions
Clinically considered surfaces that move with care.

Introducing MOVI®

A new name in pressure care.

The MOVI® range of therapeutic surfaces is designed to meet the diverse and evolving needs of patients and residents across aged care and clinical settings. Engineered for performance, safety, and infection control, each MOVI® product supports effective pressure injury prevention, improved comfort, and operational ease.

At the core of the MOVI range is the gentle alternating pressure therapy — delivering precise, low-impact pressure cycles that minimise shear, protect fragile skin, and enhance comfort —offering superior pressure care support for high-risk patients.

From overlays to powered systems, MOVI® surfaces enable a scalable approach to pressure care adapting as care needs change.

MOVI® Air T2 Overlay



Affordable clinical support.

The MOVI® Air T2 Overlay is designed to improve clinical outcomes by providing an additional layer of pressure redistribution. Lightweight and easy to deploy, it optimizes patient care for low-to-moderate risk residents.

Features

- Gentle alternating pressure therapy reduces shear, protects skin and enhances comfort
- Auto-sensing technology intelligently calibrates air pressure in real-time to the patient weight, ensuring optimal pressure redistribution and patient support
- Therapy modes: Alternating Pressure Therapy, Static Support, Max Firm, Adjustable Cycle Time, Seat Inflation Mode
- 4-way stretch, vapour permeable and fire retardant cover with fully welded seams for optimal comfort, hygiene and safety
- One-hand CPR release, positioned at the head end for fast emergency access
- Headzone with static cells for added comfort and stability
- Secures easily to any foam mattress

Code	Type	Width	Length	Height	Max User Weight
MA4100	Single	900mm	1980mm	127mm	180kg
MA4110	King Single	1050mm	1980mm	127mm	180kg

MOVI® Air T4 Replacement Mattress



A dependable support surface for consistent clinical care.

The MOVI® Air T4 Replacement Mattress delivers full-surface clinical performance, combining durable construction with therapeutic comfort for everyday use in aged care and hospital settings. Designed to support pressure injury prevention and management, it offers consistent pressure redistribution, advanced infection control properties, and an ergonomic design that enhances safety and efficiency for clinical staff during patient care and repositioning tasks.

Features

- Gentle alternating pressure therapy reduces shear, protects skin and enhances comfort
- Auto-sensing technology intelligently calibrates air pressure in real-time to the patient weight, ensuring optimal pressure redistribution and patient support
- Therapy modes: Alternating Pressure Therapy, Static Support, Max Firm, Adjustable Cycle Time, Seat Inflation Mode
- 8" Cell-on-Cell air system enhances client immersion and envelopment while delivering continuous pressure relief and support — even during power loss or transport — making it ideal for high-risk pressure injury prevention
- 4-way stretch, vapour permeable and fire retardant cover with fully welded seams for optimal comfort, hygiene and safety
- One-hand CPR release, positioned at the head end for fast emergency access
- Headzone with static cells for added comfort and stability
- Individual cell deflation allows targeted offloading, providing localised pressure relief to protect high-risk areas and support advanced pressure injury prevention and healing

Code	Type	Width	Length	Height	Max User Weight
MA4130	Single	900mm	1980mm	200mm	200kg
MA4140	King Single	1050mm	1980mm	200mm	200kg

MOVI® Comfort C3 Pumpless Mattress



Self-adjusting air technology delivering clinical-grade pressure relief without power.

The MOVI® Comfort C3 Pumpless Mattress incorporates self-adjusting air cell technology to provide dynamic pressure redistribution without the need for power, pumps, or maintenance. Ideal for facilities aiming to reduce complexity while maintaining clinical outcomes.

Features

- Automatically adapts to patient movement and weight shifts for consistent pressure relief and optimal support
- Promotes mobility, enhances comfort, and supports pressure injury prevention strategies
- Suitable for moderate to high pressure injury risk
- Multi-layer foam construction with integrated shear-reduction layers for effective pressure redistribution and enhanced skin protection
- 4-way stretch, vapour permeable and fire retardant cover with fully welded seams for optimal comfort, hygiene and safety
- Reinforced sides for safe patient transfers and edge support

Code	Type	Width	Length	Height	Max User Weight
MA4070	Single	850mm	1980mm	200mm	230kg
MA4080	King Single	1050mm	1980mm	200mm	230kg

MOVI® Comfort C4 Hybrid Mattress



Advanced pressure therapy for high-acuity care.

The MOVI Comfort C4 Hybrid Mattress provides versatile therapy options for both the prevention and treatment of pressure injuries, offering powered and non-powered modes with alternating and self-adjusting air technologies to deliver consistent, adaptive relief and advanced clinical support. This system supports patients at very high risk of pressure injuries or those with existing wounds, combining intuitive control with a quiet, low-impact footprint.

Features

- Gentle alternating pressure therapy reduces shear, protects skin and enhances comfort
- Auto-sensing technology intelligently calibrates air pressure in real-time to the patient weight, ensuring optimal pressure redistribution and patient support
- Therapy modes: Alternating Pressure Therapy, Static Support, Max Firm, Adjustable Cycle Time, Seat Inflation Mode
- Multi-layer foam and air cell construction with integrated shear-reduction layers for effective pressure redistribution and enhanced skin protection
- 4-way stretch, vapour permeable and fire retardant cover with fully welded seams for optimal comfort, hygiene and safety

Code	Type	Width	Length	Height	Max User Weight
MA4030	Single	880mm	1980mm	180mm	230kg
MA4040	King Single	1050mm	1980mm	180mm	230kg

Why Choose MOVI®




Clinically aligned

Developed to support pressure injury prevention guidelines



Scalable care

One range, four levels of support across risk profiles



Simplified Training

Consistent interface and handling features



Easy maintenance

Durable, wipe-clean covers with welded seams





Infection Control Ready

Designed to reduce cross-contamination risks



Versatile Integration

Compatible with standard and specialist beds



Find the Right Fit for Your Facility

MOVI® surfaces provide a complete solution for pressure care, allowing facilities to match the right surface to the right resident—while preparing for future changes in care needs.

For clinical evaluations, demonstrations or to speak with a specialist

1300 930 930
product@alphalifecare.com.au



MOVI® Comfort in every move.
Support at every stage.