



ook snow SC

PRODUCT CODE BD4630 | HOSPITAL BED

Features

- 5 DC Motor / 120V system auto-contour
- Battery back-up
- 4 x dual 127mm castors - central locking
- Integrated 3 stage length extension
- 4 corner roller bumpers
- CPR release handle
- Patient control lockout
- Flat deck for ease of cleaning
- Open architecture
- Backrest angle indicator
- One handed slow release head & foot end siderails
- One step Trendelenburg and one step vascular foot position
- Compliant with standard IEC 60601-2-52
- Boostless™ backrest system and backrest angle indicator
- Patient Management System
- Built in touch screen on foot end of bed
- Built in patient monitoring and 3 sensitivity level bed exit detection with alarms
- Built in bed exit side view lights
- Built in weigh scale

Controls

- Built in nurse controls on foot end
- Built in nurse controls on exterior of head end side rails
- Built in patient controls on interior of head end side rails

Specifications

Height Range

- Low Height - 254mm
- High Height - 760mm

Mattress Surface Dimensions

- Width - 889mm
- Length - 2032mm (extendable to max 2280mm)

Overall External Dimensions

- Overall Length (including bumpers) - 2240mm
- Overall Width - 1030mm (extendable to max 2490mm)

Weight Capacity

- Safe Working Load - 272kg
- Max Patient Weight - 207kg